



## Presentation

This convenient set includes the necessary tools to clean and remove cerumen from earmolds and hearing systems.

It contains a vent brush and a cleaning brush.

The cleaning brush is equipped with a magnetic to remove easily the batteries from the hearing devices, and with the nylon bristles which allow a delicate cleaning of the hearing devices.

It is also provided with a small hook on one of the extremities to scratch delicately the dirtied orifices.

## Properties

The MultiTool meets the requirements of Class I Medical Devices, according to the Regulation 2017/745/EU related to Medical Devices, while respecting environmental and safety requirements for health.

As manufacturer of Medical Devices, MG Development makes sure of the origin of the raw materials of all the materials and the coloring agents are complying with the directives and the existing standards (**REACH** for the restriction of chemical substances).

## Precautionary statements

Keep out of reach of children.