



## Presentation

The gel moisturizes and lubricates the ear canal, allows a better tolerance of the hearing aid and prevents irritations.

The gel is presented in a 0.16, 0.23 or 0.50 fl oz (5, 7 or 15ml) plastic tube.

## Properties

Contains a natural active component of Indian Frankincense, acknowledged for its potent anti-irritant effects.

The moisturizing formula of the gel makes the fitting and removing of the hearing aid, earmold or hearing protection easier.

Skin compatibility has been tested by an Institute of clinical expertise; it does not induce any skin irritations.

It meets the requirements of Class I Medical Devices, according to the Regulation 2017/745/EU related to Medical Devices, while respecting environmental and safety requirements for health.

## Using advices

Apply the cream on the hearing aid, earmold or hearing protector every morning.

**Do not put the Gel in the openings!**

## Precautionary statements

External Use.  
Do not swallow.  
Keep out of reach of children.

## Product composition

- Water
- PPG-15 Stearyl Ether
- Glycerin
- Steareth-2 & 21
- Isopropylmyristate

